

A9572

Leader in Biomolecular Solutions for Life Science



DAB2IP Mouse mAb

Catalog No.: A9572

Basic Information

Observed MW

125kDa

Calculated MW

132kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB,IP,ELISA

Cross-Reactivity

Human,Mouse

Background

DAB2IP is a Ras (MIM 190020) GTPase-activating protein (GAP) that acts as a tumor suppressor. The DAB2IP gene is inactivated by methylation in prostate and breast cancers (Yano et al., 2005 [PubMed 15386433]).

Recommended Dilutions

WB	1:500 - 1:2000
IP	1:50 - 1:100

Immunogen Information

Gene ID

153090

Swiss Prot

Q5VWQ8

Immunogen

A synthetic peptide of human DAB2IP

Synonyms

AIP1; AIP-1; AF9Q34; DIP1/2; DAB2IP

Contact



www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.