

A9687

Leader in Biomolecular Solutions for Life Science



FER Rabbit pAb

Catalog No.: A9687

Basic Information

Observed MW

110kDa

Calculated MW

95kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene is a member of the FPS/FES family of non-transmembrane receptor tyrosine kinases. It regulates cell-cell adhesion and mediates signaling from the cell surface to the cytoskeleton via growth factor receptors. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome X.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

2241

Swiss Prot

P16591

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 200-450 of human FER (NP_005237.2).

Synonyms

TYK3; PPP1R74; p94-Fer; FER

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using FER antibody (A9687) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 10s.