

A9766

Leader in Biomolecular Solutions for Life Science



PTPRZ1 Rabbit pAb

Catalog No.: A9766

Basic Information

Observed MW

Refer to figures

Calculated MW

255kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the receptor protein tyrosine phosphatase family. Expression of this gene is restricted to the central nervous system (CNS), and it may be involved in the regulation of specific developmental processes in the CNS. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

5803

Swiss Prot

P23471

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PTPRZ1 (NP_001193767.1).

Synonyms

PTPZ; HPTPZ; PTP18; PTPRZ; RPTPB; HPTPzeta; PTP-ZETA; RPTPbeta; phosphacan; R-PTP-zeta-2; PTPRZ1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.