

A9804

Leader in Biomolecular Solutions for Life Science



## MTERFD1 Rabbit pAb

Catalog No.: A9804

### Basic Information

**Observed MW**

48kDa

**Calculated MW**

48kDa

**Category**

Mouse Monoclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human, Mouse

### Background

Enables transcription cis-regulatory region binding activity. Involved in negative regulation of transcription, DNA-templated. Located in cytosol and mitochondrion.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

51001

**Swiss Prot**

Q96E29

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 69-300 of human MTERFD1 (NP\_057026.3).

**Synonyms**

CGI-12; MTERFD1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using MTERFD1 antibody (A9804) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 30s.