

A9999

Leader in Biomolecular Solutions for Life Science



## PELI3 Rabbit pAb

Catalog No.: A9999

### Basic Information

**Observed MW**

50kDa

**Calculated MW**

51kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human

### Background

The protein encoded by this gene is a scaffold protein and an intermediate signaling protein in the innate immune response pathway. The encoded protein helps transmit the immune response signal from Toll-like receptors to IRAK1/TRAF6 complexes. Several transcript variants encoding different isoforms have been found for this gene.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

246330

**Swiss Prot**

Q8N2H9

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PELI3 (NP\_659502.2).

**Synonyms**

PELI3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data

