

Rabbit anti Myc-Tag pAb

Catalog No	AE009	Category	Tag Antibodies
Applications	WB	Observed MW	Refer to Figures
Cross-Reactivity		Calculated MW	

Immunogen Information

Immunogen Recombinant protein of MYC-tag

Gene ID

Swiss Prot

Synonyms

Product information

Source Rabbit

Isotype IgG


Purity Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:1000 - 1:10000



Western blot analysis of over-expressed Myc-tagged protein in 293T cell using Myc-tag antibody (AE009) at different dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.