

Rabbit anti Myc-Tag pAb

Catalog No.	AE009	Category	Tag Antibodies
Applications	WB	Observed MW	Refer to Figures
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen Recombinant protein of MYC-tag

Gene ID

Swiss prot

Synonyms

Product information

Source Rabbit

Isotype IgG


Purification method Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:1000 - 1:10000



Western blot - Myc-tag antibody (AE009)