

## Rabbit anti Myc-Tag pAb

|                         |       |                      |                  |
|-------------------------|-------|----------------------|------------------|
| <b>Catalog No.</b>      | AE009 | <b>Category</b>      | Tag Antibodies   |
| <b>Applications</b>     | WB    | <b>Observed MW</b>   | Refer to Figures |
| <b>Cross-reactivity</b> |       | <b>Calculated MW</b> |                  |

### Immunogen Information

|                   |                                |
|-------------------|--------------------------------|
| <b>Immunogen</b>  | Recombinant protein of MYC-tag |
| <b>Gene ID</b>    |                                |
| <b>Swiss prot</b> |                                |
| <b>Synonyms</b>   |                                |


### Product information

|                            |  |
|----------------------------|--|
| <b>Source</b>              | Rabbit   |
| <b>Isotype</b>             | IgG  |
| <b>Purification method</b> | Affinity purification  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |

### Background

### Recommended Dilutions

WB 1:1000 -  
1:10000



Western blot - Rabbit anti Myc-Tag pAb (AE009)