

Rabbit anti GFP-Tag pAb

Catalog No	AE011	Category	Tag Antibodies
Applications	WB,IHC,IF	Observed MW	
Cross-Reactivity		Calculated MW	

Immunogen Information

Immunogen	Recombinant protein of GFP
Gene ID	
Swiss Prot	
Synonyms	


Product information

Source	Rabbit
Isotype	IgG
Purity	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB	1:2000 - 1:3000
IHC	1:50 - 1:200
IF	1:50 - 1:200



Western blot analysis of over-expressed GFP fusion protein in 293 cell using GFP-Tag antibody (AE011) at different dilution. Each lane was loaded with 2 ug cell lysate. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.