

Rabbit anti GFP-Tag pAb

Catalog No.	AE011	Category	Tag Antibodies
Applications	WB, IHC, IF	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Recombinant protein of GFP
Gene ID	
Swiss prot	
Synonyms	

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

Recommended Dilutions


WB 1:2000 -
1:3000

IHC 1:50 -
1:200

IF 1:50 -
1:200



Western blot - Rabbit anti GFP-Tag pAb (AE011)



Immunofluorescence - Rabbit anti GFP-Tag pAb (AE011)