

Rabbit anti GFP-Tag pAb

Catalog No.	AE011	Category	Tag Antibodies
Applications	WB, IHC, IF	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen Recombinant protein of GFP

Gene ID

Swiss prot

Synonyms

Product information

Source Rabbit

Isotype IgG

Purification method Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

Recommended Dilutions

WB 1:2000 - 1:3000

IHC 1:50 - 1:200

IF 1:50 - 1:200



Western blot analysis of over-expressed GFP fusion protein in 293 cell using GFP-Tag antibody (AE011) at different dilution. Each lane was loaded with 2 ug cell lysate.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.