

AE087

Leader in Biomolecular Solutions for Life Science



Magnetic Beads-conjugated Mouse anti HA-Tag mAb

Catalog No.: AE087

Basic Information

Observed MW

52kDa

Calculated MW

Category

SMab Recombinant Monoclonal Antibody

Applications

IP, ELISA

Cross-Reactivity

Species independent

Background

Recommended Dilutions

IP 1:50 - 1:200

Immunogen Information

Gene ID

Swiss Prot

Immunogen

A synthetic peptide corresponding to HA tag.

Synonyms

Contact



www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.