

AP0510

Leader in Biomolecular Solutions for Life Science



Phospho-Fyn-S21 Rabbit pAb

Catalog No.: AP0510

2 Publications

Basic Information

Observed MW

60kDa

Calculated MW

61kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

2534

Swiss Prot

P06241

Immunogen

A synthetic phosphorylated peptide around S21 of human Phospho-Fyn-S21 (NP_002028.1).

Synonyms

SLK; SYN; p59-FYN; Phospho-Fyn-S21

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of lysates from A-431 cells, using Phospho-Fyn-S21 Rabbit pAb (AP0510) at 1:1000 dilution. HL-60 cells were treated by 10% FBS after serum-starvation overnight.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% BSA.