

AP1057

Leader in Biomolecular Solutions for Life Science



Phospho-Paxillin-Y118 Rabbit pAb

Catalog No.: AP1057

Basic Information

Observed MW

Refer to figures

Calculated MW

65kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5829

Swiss Prot

P49023

Immunogen

A synthetic phosphorylated peptide around Y118 of human Paxillin (NP_001074324.1?).

Synonyms

Paxillin; PXN; paxillin; Phospho-Paxillin-Y118

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.