

AP1108

Leader in Biomolecular Solutions for Life Science



Phospho-STK11-S428 Rabbit pAb

Catalog No.: AP1108

Basic Information

Observed MW

65kDa

Calculated MW

49kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene is a serine/threonine kinase that regulates cell polarity and energy metabolism and functions as a tumor suppressor. Mutations in this gene have been associated with the autosomal dominant Peutz-Jeghers syndrome, as well as with skin, pancreatic, and testicular cancers.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

6794

Swiss Prot

Q15831

Immunogen

A synthetic phosphorylated peptide around S428 of human STK11 (NP_000446.1).

Synonyms

PJS; LKB1; hLKB1; Phospho-STK11-S428

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of lysates from NIH/3T3 cells, using Phospho-STK11-S428 Rabbit pAb (AP1108) at 1:1000 dilution or LKB1 antibody (A2122). NIH/3T3 cells were treated by PMA/TPA (200 nM) at 37°C for 30 minutes after serum-starvation overnight.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% BSA.

Detection: ECL Basic Kit (RM00020).

Exposure time: 60s.