

AP1128

Leader in Biomolecular Solutions for Life Science



## Phospho-DAPP1-Y139 Rabbit pAb

Catalog No.: AP1128

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

32kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse

### Background

Enables phosphatidylinositol-3,4,5-trisphosphate binding activity and phosphatidylinositol-3,4-bisphosphate binding activity. Predicted to be involved in signal transduction. Located in plasma membrane.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

27071

#### Swiss Prot

Q9UN19

#### Immunogen

A synthetic phosphorylated peptide around Y139 of human DAPP1 (NP\_055210.2).

#### Synonyms

BAM32; Phospho-DAPP1-Y139

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.