

AP1203

Leader in Biomolecular Solutions for Life Science



Phospho-VANGL2-S11/S14 Rabbit mAb

Catalog No.: AP1203

Recombinant

Basic Information

Observed MW

Refer to figures

Calculated MW

60kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

CloneNo number

ARC5007-02

Background

The protein encoded by this gene is a membrane protein involved in the regulation of planar cell polarity, especially in the stereociliary bundles of the cochlea. The encoded protein transmits directional signals to individual cells or groups of cells in epithelial sheets. This protein is also involved in the development of the neural plate.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

57216

Swiss Prot

Q9ULK5

Immunogen

A phospho specific peptide corresponding to residues surrounding S11/S14 of human VANGL2

Synonyms

LPP1; LTAP; STB1; STBM; STBM1; Phospho-VANGL2-S11/S14

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 0.05% BSA, 50% glycerol, pH7.3.