

AP1205

Leader in Biomolecular Solutions for Life Science



## Phospho-VANGL2-T78/S79 Rabbit mAb

Catalog No.: AP1205

Recombinant

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

60kDa

#### Category

SMab Recombinant Monoclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

#### CloneNo number

ARC5015-09

### Background

The protein encoded by this gene is a membrane protein involved in the regulation of planar cell polarity, especially in the stereociliary bundles of the cochlea. The encoded protein transmits directional signals to individual cells or groups of cells in epithelial sheets. This protein is also involved in the development of the neural plate.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

57216

#### Swiss Prot

Q9ULK5

#### Immunogen

A phospho specific peptide corresponding to residues surrounding T78/S79 of human VANGL2

#### Synonyms

LPP1; LTAP; STB1; STBM; STBM1; Phospho-VANGL2-T78/S79

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 0.05% BSA, 50% glycerol, pH 7.3.