

AP1379

Leader in Biomolecular Solutions for Life Science



# Phospho-FOXO1-T24 + FOXO3A-T32 + FOXO4-T28 Rabbit mAb

Catalog No.: AP1379

Recombinant

## Basic Information

### Observed MW

78-82kDa

### Calculated MW

70kDa/71kDa/48kDa/54kDa

### Category

SMab Recombinant Monoclonal Antibody

### Applications

WB, ELISA

### Cross-Reactivity

Human, Mouse

### CloneNo number

ARC57268

## Background

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma. [provided by RefSeq, Jul 2008]

## Recommended Dilutions

WB 1:1000 - 1:5000

## Immunogen Information

### Gene ID

2308/2309/4303

### Swiss Prot

Q12778/O43524/P98177

### Immunogen

A synthetic phosphorylated peptide around T24 & T32 & T28 of human FOXO1/FOXO3A/FOXO4FOXO1.

### Synonyms

FoxO1/FoxO3a/FoxO4; Phospho-FOXO1-T24 + FOXO3A-T32 + FOXO4-T28

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

## Validation Data



Western blot analysis of Jurkat, using Phospho-FOXO1-T24 + FOXO3A-T32 + FOXO4-T28 antibody (AP1379) at 1:2000 dilution. Jurkat cells were treated by Calyculin A (100 nM) at 37°C for 30 minutes. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 20s.



Western blot analysis of NIH/3T3, using Phospho-FOXO1-T24 + FOXO3A-T32 + FOXO4-T28 antibody (AP1379) at 1:2000 dilution. NIH/3T3 cells were treated by Calyculin A (100 nM) at 37°C for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 20s.