

FITC Goat Anti-Mouse IgG (H+L)

Catalog No.	AS001	Category	Secondary Antibodies
Applications	IF, FC	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Mouse IgG
Gene ID	
Swiss prot	
Synonyms	

Product information

Source	Goat
Isotype	Fluorescein conjugated IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

ICC	1:20 - 1:100
FC	1:20 - 1:100

