

FITC Goat Anti-Mouse IgG (H+L)

Catalog No.	AS001	Category	Secondary Antibodies
Applications	IF, FC	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen Mouse IgG

Gene ID

Swiss prot

Synonyms

Product information

Source Goat

Isotype Fluorescein conjugated IgG

Purification method Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

ICC 1:20 - 1:100

FC 1:20 - 1:100