

HRP Goat Anti-Human IgG (H+L)

Catalog No.	AS002	Category	Secondary Antibodies
Applications	WB, IHC	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Human IgG
Gene ID	
Swiss prot	
Synonyms	

Product information

Source	Goat
Isotype	Horseradish peroxidase conjugated IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB	1:2000 - 1:50000
IHC	1:200 - 1:2000
ICC	1:200 - 1:2000