

HRP Goat Anti-Mouse IgG (H+L)

Catalog No	AS003	Category	Secondary Antibodies
Applications	Elisa, WB	Observed MW	
Cross-Reactivity		Calculated MW	

Immunogen Information

Immunogen	Mouse IgG
Gene ID	
Swiss Prot	
Synonyms	

Product information

Source	Goat
Isotype	Horseradish peroxidase conjugated IgG
Purity	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

Elisa	1:5000 - 1:10000
WB	1:2000 - 1:10000