

HRP Goat Anti-Mouse IgG (H+L)

Catalog No.	AS003	Category	Secondary Antibodies
Applications	Elisa, WB	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Mouse IgG
Gene ID	
Swiss prot	
Synonyms	

Product information

Source	Goat
Isotype	Horseradish peroxidase conjugated IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: 50% glycerol, pH7.3.

Background

Recommended Dilutions

Elisa 1:5000 -
1:10000
WB 1:2000 -
1:10000

