

AMCA Goat Anti-Rabbit IgG (H+L)

Catalog No.	AS020	Category	Secondary Antibodies
Applications	ICC, FC	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen Rabbit IgG

Gene ID

Swiss prot

Synonyms

Product information

Source Goat

Isotype AMCA conjugated IgG

Purification method Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

ICC 1:50 - 1:200

FC 1:50 - 1:200