

AMCA Goat Anti-Rabbit IgG (H+L)

| | | | |
|-------------------------|--------|----------------------|----------------------|
| Catalog No | AS020 | Category | Secondary Antibodies |
| Applications | ICC,FC | Observed MW | |
| Cross-Reactivity | | Calculated MW | |

Immunogen Information

Immunogen Rabbit IgG

Gene ID

Swiss Prot

Synonyms

Product information

Source Goat

Isotype AMCA conjugated IgG

Purity Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

ICC 1:50 - 1:200

FC 1:50 - 1:200