

## HRP Goat Anti-Guinea Pig IgG (H+L)

<b>Catalog No.</b>	AS025	<b>Category</b>	Secondary Antibodies
<b>Applications</b>	WB, IHC	<b>Observed MW</b>	
<b>Cross-reactivity</b>		<b>Calculated MW</b>	

### Immunogen Information

<b>Immunogen</b>	Guinea Pig IgG
<b>Gene ID</b>	
<b>Swiss prot</b>	
<b>Synonyms</b>	

### Product information

<b>Source</b>	Goat
<b>Isotype</b>	Horseradish peroxidase conjugated IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

### Recommended Dilutions

WB	1:2000 - 1:50000
IHC	1:200 - 1:2000
ICC	1:200 - 1:2000

