

HRP Rabbit Anti-Goat IgG (H+L)

Catalog No.	AS029	Category	Secondary Antibodies
Applications	WB, IHC	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Goat IgG
Gene ID	
Swiss prot	
Synonyms	

Product information

Source	Rabbit
Isotype	Horseradish peroxidase conjugated IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:2000 -
1:50000

IHC 1:200 -
1:2000

ICC 1:200 -
1:2000

