

HRP Rabbit Anti-Goat IgG (H+L)

| | | | |
|-------------------------|---------|----------------------|----------------------|
| Catalog No. | AS029 | Category | Secondary Antibodies |
| Applications | WB, IHC | Observed MW | |
| Cross-reactivity | | Calculated MW | |

Immunogen Information

| | |
|-------------------|----------|
| Immunogen | Goat IgG |
| Gene ID | |
| Swiss prot | |
| Synonyms | |

Product information

| | |
|----------------------------|--|
| Source | Rabbit |
| Isotype | Horseradish peroxidase conjugated IgG |
| Purification method | Affinity purification |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |

Background

Recommended Dilutions

| | |
|-----|------------------|
| WB | 1:2000 - 1:50000 |
| IHC | 1:200 - 1:2000 |
| ICC | 1:200 - 1:2000 |

