

Monkey IFN-gamma ELISA Kit

Catalog NO.: RK00161

version: 2.0

This package insert must be read in its entirety before using this product

1 www.abclonal.com



Introduction

The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of IFN-gamma in Monkey serum, plasma, cell culture supernatants and other biological fluids.

Principle of the Assay

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for Monkey IFN-gamma has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IFN-gamma present is bound by the immobilized antibody. After washing away any unbound substances, a biotin-conjugated detection antibody specific for Monkey IFN-gamma is added to the wells. After incubating, an enzyme-linked streptavidin is added and binds to biotin. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IFN-gamma bound in the initial step. The color development is stopped, and the absorbance is measured.



Material Provided & Storage Conditions

Unopened kits can be stored at $2-8^{\circ}$ C for 1 year, and opened products must be used within 1 month.

Part	Size	Cat. No.	Storage of opened/reconstituted material
Antibody Coated Plate	8×12	RM00672	Put the unused slats back in the aluminum foil bag with the desiccant and reseal them. They can be stored at 2-8° C for 1 month.
Standard 2 vials		RM00669	It is not recommended to use again after redissolving.
Concentrated Biotin Conjugate Antibody (30×)	1 ×400ul	RM00670	Store at 2-8° c for 1
Streptavidin-HRP Concentrated (30 ×)	1 ×400u1	RM00671	Store at 2-8° c for 1



Standard/Sample	1 ×20mL	RM00023	
Diluent (R1)			
Biotin-Conjugate	4 1/40 7	Divosos (
Antibody Diluent	1 ×12mL	RM00024	Store at 2-8° c for 1
(R2)			month *
Streptavidin-HRP	1 × 12mL	RM00025	month *
Diluent(R3)	1 /(12mL		
Wash Buffer(20x)	1 ×30mL	RM00026	
TMB Substrate	1 ×12mL	RM00027	
Stop Solution	1 ×6mL	RM00028	
Plate Sealers	4 Strips		
Specification	1		



Other Supplies Required

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 630 nm or 570 nm.
- 2. Pipettes and pipette tips.
- 3. Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- Incubator.
- 6. Test tubes for dilution of standards and samples.



Precautions

* FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

- Any variation in diluent, operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.
- Variations in sample collection, processing, and storage may cause sample value differences.
- Reagents may be harmful, if ingested, rinse it with an excess amount of tap water.
- Stop Solution contains strong acid. Wear eye, hand, and face protection.
- Please perform simple centrifugation to collect the liquid before use.
- Do not mix or substitute reagents with those from other lots or other sources.
- Adequate mixing is particularly important for good result.
 Use a mini-vortexer at the lowest frequency.
- Mix the sample and all components in the kits adequately, and use clean plastic container to prepare all diluents.
- 9. Both the sample and standard should be assayed in



duplicate, and reagents should be added in sequence in accordance with the requirement of the specification.

- 10. Reuse of dissolved standard is not recommended.
- 11. The kit should not be used beyond the expiration date on the kit label.
- The kit should be away from light when it is stored or incubated.
- 13. To reduce the likelihood of blood-borne transmission of infectious agents, handle all serum, plasma, and other biological fluids in accordance with NCCLS regulations.
- To avoid cross contamination, please use disposable pipette tips.
- 15. Please prepare all the kit components according to the Specification. If the kits will be used several times, please seal the rest strips and preserve with desiccants. Do use up within 2 months.
- 16. This assay is designed to eliminate interference by other factors present in biological samples.
- 17. Until all factors have been tested in this assay, the possibility of interference cannot be excluded.
- 18. The 48T kit is also suitable for the specification.



Sample Collection & Storage

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernatant: Remove particulates by centrifugation. Assay immediately or aliquot and store samples at \leq -20 °C. Avoid repeated freeze-thaw cycles.

Serum: Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 xg. Remove serum and assay immediately or aliquot and store samples at \leq -20 °C. Avoid repeated freeze-thaw cycles.

Plasma: Collect plasma using EDTA or Heparin as an anticoagulant. Centrifuge for 15 minutes at $1000 \times g$ within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. (Note: Citrate plasma has not been validated for use in this assay.)

Other biological fluids: Centrifuge samples for 20 minutes at $1,000\times g$. Collect the supernatants and assay immediately or store samples in aliquot at -20° C or -80° C for later use. Avoid



repeated freeze-thaw cycles.

Note: It is suggested that all samples in one experiment be collected at the same time of the day. Avoid hemolytic and hyperlipidemia sample for serum and plasma.

Reagent Preparation


Bring all reagents to room temperature before use. If crystals have formed in the concentrate, Bring the reagent to room temperature and mix gently until the crystals have completely dissolved.

Standard - Reconstitute the Standard Lyophilized with 1.0 mL Standard/Sample Diluent(R1). This reconstitution produces a stock solution of 500pg/mL. Mix the standard to ensure complete reconstitution and allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Use the 500pg/mL standard stock to produce a dilution series

9 www.abclonal.com



(below) with Standard/Sample Diluent(R1). Mix each tube thoroughly and change pipette tips between each transfer (recommended concentration for standard curve: 250, 125, 62.5, 31.25, 15.6, 7.8, 3.9, 0pg/mL). Use diluted standards within 60 minutes of preparation.



Working Biotin Conjugate Antibody - Dilute 1:30 of Concentrated Biotin Conjugate Antibody (30x) with Biotin-Conjugate Antibody Diluent (R2) before use. For example: Add 66 µL of Concentrated Biotin Conjugate Antibody (30x) to 1914 µL Biotin-Conjugate Antibody Diluent (R2) solution for the preparation of 1980 µL Working Biotin-Conjugate Antibody Buffer.

Working Streptavidin-HRP - Dilute 1:30 of Concentrated Streptavidin-HRP (30x) with Streptavidin-HRP Diluent (R3)



before use. For example: Add 66 μ L of Concentrated Streptavidin-HRP (30x) to 1914 μ L Streptavidin-HRP Diluent (R3) to prepare 1980 μ L Working Streptavidin-HRP Buffer.

Wash Buffer - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 1:20 with double distilled or deionized water before use. For example: Add 20 mL of Wash Buffer Concentrate to 380 mL of deionized or distilled water to prepare 400 mL of Wash Buffer.

Assay Procedure

11 www.abclonal.com



Bring all reagents and samples to room temperature before use. It is recommended that all standards, controls, and samples be assayed in duplicate.

- Prepare all reagents, working standards, and samples as directed in the previous sections.
- Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
- 3. Add wash buffer 350 μ L/well, aspirate each well after holding 40 seconds, repeating the process two times for a total of three washes.
- 4. Add 100 μL Standard/sample Diluent (R1) in a blank well.
- 5. Add 100 µL different concentration of standard or sample in other wells, Cover with the adhesive sealer provided. Incubate for 2 hours at 37°C. Record the plate layout of standards and sample assay.
- Prepare the Concentrated Biotin Conjugate Antibody (30x)
 Working Solution 15 minutes early before use.
- 7. Repeat the aspiration/wash as in step 3.
- Add 100 μL Working Biotin Conjugate Antibody in each well, cover with new adhesive sealer provided. Incubate for 1 hour at 37°C.



- Prepare the Streptavidin-HRP Concentrated (30x) Working Solution 15minutes early before use.
- 10. Repeat the aspiration/wash as in step 3.
- Add 100 µL Working Streptavidin-HRP in each well, cover with new adhesive sealer provided. Incubate for 0.5 hour at 37°C.
- 12. Repeat the aspiration/wash as in step 3.
- During the incubation, turn on the microplate reader to warm up for 30 minutes before measuring.
- Add 100 μL TMB Substrate to each well. Incubate for 15-20 minutes at 37°C. Protect from light.
- 15. Add 50 μL Stop Solution, determine the optical density of each well within 5 minutes, using a Microplate reader set to 450 nm. If wavelength correction is available, set to 570 nm or 630 nm. If wavelength correction is not available, subtract readings at 570 nm or 630 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may cause higher value and less accurate result.

Assay Procedure Summary



Prepare the standard and reagents

Wash 3 times

Ţ

Add 100ul of standards or test samples to each well Incubate for 2 hours at 37°C, then wash 3 times

1

Add 100ul Working Biotin Conjugate Antibody Incubate for 1 hours at $37\,^{\circ}\mathrm{C}$, then wash 3 times

ļ

Add 100ul Working Streptavidin-HRP

Incubate for 0.5 hour at 37°C , then wash 3 times

1

Add 100ul Substrate Solution

Incubate for 15-20 min at 37°C under dark condition

ļ

Add 50ul Stop Solution

ı

Detect the optical density within 5 minutes under 450nm.

Correction Wavelength set at 570nm or 630nm

Calculation of Results


14 www.abclonal.com



- Average the duplicate readings for each standard, control and sample, and subtract the average zero standard optical density (0.D.).
- 2. Create a standard curve by reducing the data using computer software capable of generating a four- parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the Y-axis against the concentration on the X-axis and draw a best fit curve through the points on a log/log graph. The data may be linearized by plotting the log of the IFN-gamma concentrations versus the log of the 0.D. on a linear scale, and the best fit line can be determined by regression analysis.
- If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.



Typical Data



The standard curves are provided for demonstration only. A standard curve should be generated for each set of IFN-gamma assayed.

Detection Range

 $3.\,9\text{--}250 \mathrm{pg/mL}$



Sensitivity

The minimum detectable dose (MDD) of IFN-gamma typically less than 1.95pg/mL. The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration

Specificity

This assay recognizes both recombinant and natural Monkey IFN-gamma. The factors listed below were prepared at 50ng/ml and assayed for cross-reactivity. No significant cross-reactivity was observed with the following:

Recombinant Monkey: Recombinant Pig:

IL-1 IFN- γ

IL-2

IL-3

IL-4

IL-5

IL-6



IL-10

II.-12

TNF

SCF

LIF

GM-CSF

Note:

Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between IFNgamma and all the analogues, therefore, cross reaction may still exist.



Precision

Intra-plate Precision

The three samples with known concentrations were repeated 20 times on the same board to calculate the coefficient of variation (CV) for the concentrations.

Intra-Assay: CV<10%

Inter-plate Precision

The coefficient of variation (CV) of three known concentrations was calculated by 20 tests of two kits on each of the three samples.

Inter-Assay: CV<15%

	Intra-Assay Precision			Inter-	Assay Pre	cision
Sample	1	2	3	1	2	3
n	20	20	20	20	20	20
Mean(pg/mL)	56	110	240	60	120	250
Standard deviation	2. 52	5. 06	11. 76	3. 48	7. 68	14. 75
CV (%)	4. 5	4.6	4. 9	5.8	6. 4	5. 9



Recovery

Matrices listed below were spiked with certain level of IFN-gamma and the recovery rates were calculated by comparing the measured value to the expected amount of IFN-gamma in samples.

Sample	Average Recovery (%)	Range (%)	
Cell Culture Media(n=5)	97	83-103	
Serum(n=5)	104	91–106	



Linearity

The linearity of the kit was assayed by testing samples spiked with appropriate concentration of IFN-gamma and their serial dilutions. The results were demonstrated by the percentage of calculated concentration to the expected.

/	/	Cell Culture Media(n=5)	Serum(n=5)
1:2	Average of Expected (%)	100	95
	Range (%)	85-112	80-98
1:4	Average of Expected (%)	101	92
	Range (%)	92-104	87-97
1:8	Average of Expected (%)	95	88
	Range (%)	92-99	83-92
1:16	Average of Expected (%)	98	90
	Range (%)	90-106	84-96

Trouble Shooting



Problem	Possible Cause	Solution
High Background	Insufficient washing	Sufficiently wash plates as required. Ensure appropriate duration and number of washes. Ensure appropriate volume of wash buffer in each well.
	Incorrect incubation procedure	Check whether the duration and temperature of incubation are set up as required.
	Cross-contamination of samples and reagents	Be careful of the operations that could cause cross-contamination. Use fresh reagents and repeat the tests.
No signal or weak signal	Incorrect use of reagents	Check the concentration and dilution ratio of reagents. Make sure to use reagents in proper order.
	Incorrect use of	Warm the reader up before use. Make sure to set up appropriate main wavelength and correction wavelength.
	Insufficient colour reaction time	Optimum duration of colour reaction should be limited to 15-25 minutes.
	Read too late after stopping the colour reaction	Read the plate in 5 minutes after stopping the reaction.
	Matrix effect of samples	Use positive control.



	Contamination of TMB	Check if TMB substrate solution turns
Too much	substrate	blue. Use new TMB substrate solution.
	Plate sealers reused	Use a fresh new sealer in each step of experiments.
	Protein concentration in sample is too high	Do pre-test and dilute samples in optimum dilution ratio.
Poor Duplicates	Uneven addition of samples	Check the pipette. Periodically calibrate the pipette.
	Impurities and precipitates in samples	Centrifuge samples before use.
	Inadequate mixing of reagents	Mix all samples and reagents well before loading.

^{*}For research purposes only. Not for therapeutic or diagnostic purposes.