

Phospho-ERK1-T202/Y204 + ERK2-T185/Y187 Antibody Kit

Catalog No.	RK04105	Category	Antibody Kit
Applications	WB, IP	Observed MW	42kDa/44kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	36kDa/41kDa/38kDa/40kDa/43kDa

Immunogen Information

Immunogen A synthetic phosphorylated peptide around T185/Y187 of human ERK2 (NP_002736.3).
Recombinant Protein of human MAPK1/MAPK3

Gene ID 5594/5595

Swiss prot P28482/P27361

Synonyms

Product information

Source Rabbit


Isotype IgG

Purification method Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions



Western blot - Phospho-ERK1-T202/Y204 + ERK2-T185/Y187 pAb (AP0472)