

RK05648

Leader in Biomolecular Solutions for Life Science



[One Step] GDA Antibody Kit

Catalog No.: RK05648

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

51kDa

Calculated MW

51kDa

Category

Antibody kit

Background

This gene encodes an enzyme responsible for the hydrolytic deamination of guanine. Studies in rat ortholog suggest this gene plays a role in microtubule assembly. Multiple transcript variants encoding different isoforms have been found for this gene.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05648-1	GDA Rabbit pAb	51kDa
RK05648-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	51kDa

Immunogen Information

Gene ID

9615

Swiss Prot

Q9Y2T3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 165-454 of human GDA (NP_004284.1).

Synonyms

GAH; CYPIN; GUANASE; NEDASIN; GDA