

RK05652

Leader in Biomolecular Solutions for Life Science



[One Step] HACE1 Antibody Kit

Catalog No.: RK05652

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

102kDa

Calculated MW

102kDa

Category

Antibody kit

Background

This gene encodes a HECT domain and ankyrin repeat-containing ubiquitin ligase. The encoded protein is involved in specific tagging of target proteins, leading to their subcellular localization or proteasomal degradation. The protein is a potential tumor suppressor and is involved in the pathophysiology of several tumors, including Wilm's tumor.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05652-1	HACE1 Rabbit pAb	102kDa
RK05652-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	102kDa

Immunogen Information

Gene ID

57531

Swiss Prot

Q8IYU2

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HACE1 (NP_065822.2).

Synonyms

SPPRS; HACE1