

RK05670

Leader in Biomolecular Solutions for Life Science



[One Step] NDUFB10 Antibody Kit

Catalog No.: RK05670

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

21kDa

Calculated MW

21kDa

Category

Antibody kit

Background

Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrial inner membrane. Part of mitochondrial respiratory chain complex I. Implicated in nuclear type mitochondrial complex I deficiency 35.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05670-1	NDUFB10 Rabbit pAb	21kDa
RK05670-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	21kDa

Immunogen Information

Gene ID

4716

Swiss Prot

O96000

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human NDUFB10 (NP_004539.1).

Synonyms

PDSW; MC1DN35; NDUFB10