

RK05679

Leader in Biomolecular Solutions for Life Science



[One Step] PHB Antibody Kit

Catalog No.: RK05679

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

30kDa

Calculated MW

30kDa

Category

Antibody kit

Background

This gene is evolutionarily conserved, and its product is proposed to play a role in human cellular senescence and tumor suppression. Antiproliferative activity is reported to be localized to the 3' UTR, which is proposed to function as a trans-acting regulatory RNA. Several pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05679-1	Prohibitin Rabbit pAb	30kDa
RK05679-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	30kDa

Immunogen Information

Gene ID

5245

Swiss Prot

P35232

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 73-272 of human Prohibitin (NP_002625.1).

Synonyms

PHB; HEL-215; HEL-S-54e; Prohibitin