

RK05681

Leader in Biomolecular Solutions for Life Science



[One Step] POLR2H Antibody Kit

Catalog No.: RK05681

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

17kDa

Calculated MW

17kDa

Category

Antibody kit

Background

The three eukaryotic RNA polymerases are complex multisubunit enzymes that play a central role in the transcription of nuclear genes. This gene encodes an essential and highly conserved subunit of RNA polymerase II that is shared by the other two eukaryotic DNA-directed RNA polymerases, I and III. Alternative splicing results in multiple transcript variants of this gene.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05681-1	POLR2H Rabbit pAb	17kDa
RK05681-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	17kDa

Immunogen Information

Gene ID

5437

Swiss Prot

P52434

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human POLR2H (NP_006223.2).

Synonyms

RPB8; RPB17; RPABC3; POLR2H