

RK05687

Leader in Biomolecular Solutions for Life Science



[One Step] PSMA4 Antibody Kit

Catalog No.: RK05687

1 Publications

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

28kDa

Calculated MW

29kDa

Category

Antibody kit

Background

This gene encodes a core alpha subunit of the 20S proteasome, which is a highly ordered ring-shaped structure composed of four rings of 28 non-identical subunits. Proteasomes cleave peptides in an ATP- and ubiquitin-dependent manner.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05687-1	PSMA4 Rabbit pAb	29kDa
RK05687-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	28kDa

Immunogen Information

Gene ID

5685

Swiss Prot

P25789

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human PSMA4 (NP_001096137.1).

Synonyms

HC9; PSC9; HsT17706; PSMA4