

RK05734

Leader in Biomolecular Solutions for Life Science



# [One Step] UBE2I antibody Kit

Catalog No.: RK05734

## Basic Information

### Applications

### Cross-Reactivity

Human, Mouse, Rat

### Observed MW

18kDa

### Calculated MW

18kDa

### Category

Antibody kit

## Background

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Four alternatively spliced transcript variants encoding the same protein have been found for this gene.

## Product Information

### Source

Rabbit

### Purification

### Storage

Avoid repeated freeze-thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05734-1	UBE2I Rabbit pAb	18kDa
RK05734-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	18kDa

## Immunogen Information

### Gene ID

7329

### Swiss Prot

P63279

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human UBE2I (NP\_003336.1).

### Synonyms

P18; UBC9; C358B7.1; UBE2I