

RK05738

Leader in Biomolecular Solutions for Life Science



[One Step] PPP4C antibody Kit

Catalog No.: RK05738

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

36kDa

Calculated MW

35kDa

Category

Antibody kit

Background

Enables protein serine/threonine phosphatase activity. Involved in regulation of double-strand break repair via homologous recombination. Located in cytosol; nucleoplasm; and plasma membrane.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05738-1	PPP4C Rabbit pAb	35kDa
RK05738-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	36kDa

Immunogen Information

Gene ID

5531

Swiss Prot

P60510

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-307 of human PPP4C (NP_002711.1).

Synonyms

PP4; PPX; PP-X; PP4C; PPH3; PPP4; PPP4C