

RP00021

Leader in Biomolecular Solutions for Life Science



Recombinant Human BID Protein

Catalog No.: RP00021

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	637	P55957

Tags

C-His

Synonyms

BID;FP497

Product Information

Source

<I>E. coli</I>

Purification

> 95% by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of 50mM Tris, 150mM NaCl, pH 8.0. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Contact



www.abclonal.com

Background

The BH3 interacting domain death agonist (BID) is a pro-apoptotic member of the Bcl-2 protein family, which contains only the BH3 domain, and is required for its interaction with the Bcl-2 family proteins and for its pro-death activity. BID is important to cell death mediated by these proteases and thus is the sentinel to protease-mediated death signals. It is a mediator of mitochondrial damage induced by caspase-8 (CASP8); CASP8 cleaves this encoded protein, and the COOH-terminal part translocates to mitochondria where it triggers cytochrome c release.

Basic Information

Description

Recombinant Human BID Protein is produced by <I>E. coli</I> expression system. The target protein is expressed with sequence (Met1-Asp195) of human BID (Accession #NP_001187.1) fused with a 6×His tag at the C-terminus.


Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human BID at 2 μg/mL (100 μL/well) can bind BCL-XL with a linear range of 0.195-4.94 ng/mL.


Storage

Store the lyophilized protein at -20°C to -80 °C for long term.
After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human BID Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 22 kDa.



Immobilized Human BID at 2 $\mu\text{g/mL}$ (100 μL /well) can bind BCL-XL with a linear range of 0.195-4.94 ng/mL.