

RP00210

Leader in Biomolecular Solutions for Life Science



Recombinant Human EGFR Protein

Catalog No.: RP00210

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	1956	P00533

Tags

C-His

Synonyms

EGFR; ERBB; ERBB1; HER1; NISBD2; PIG61; mENA; epidermal growth factor receptor; ERBB; ERBB1; HER1; NISBD2; PIG61; mENA

Product Information

Source	Purification
HEK293 cells	> 92% by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Background

The EGFR (epidermal growth factor receptor) subfamily of receptor tyrosine kinases comprises four members: EGFR (also known as HER-1, ErbB1, or ErbB), ErbB2 (Neu, HER-2), ErbB3 (HER-3), and ErbB4 (HER-4). EGFR protein is type I transmembrane glycoprotein that binds a subset of EGF family ligands including EGF, amphiregulin, TGF-α, betacellulin, etc. The human EGFR cDNA encodes a 1210 amino acid (aa) precursor with a 24 aa signal peptide, a 621 aa extracellular domain (ECD), a 23 aa transmembrane segment, and a 542 aa cytoplasmic domain. EGFR signaling has been shown to exert action on carcinogenesis and disease progression, and thus EGFR protein is proposed as a target for cancer therapy currently, which is overexpressed in a wide variety of tumors and is the target of several anti-cancer drugs.

Basic Information

Description

Recombinant Human EGFR Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Ser645) of human EGFR (Accession #NP_005219.2.) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human EGFR at 5 μg/mL (100 μL/well) can bind Human EGF with a linear range of 7-25 ng/mL.

Storage


Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact




www.abclonal.com

Validation Data



Recombinant Human EGFR Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 110-115kDa.



Immobilized Human EGFR at 5 $\mu\text{g/mL}$ (100 $\mu\text{L/well}$) can bind Human EGF with a linear range of 7-25ng/mL.