

RP00250

Leader in Biomolecular Solutions for Life Science



# Recombinant Human Ephrin-B1/EFNB1 Protein

Catalog No.: RP00250

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	1947	P98172

### Tags

C-hFc&His

### Synonyms

EFNB1;CFND;CFNS;EFB1;EFL3;EPLG2;Elk-L;LERK2;ephrin-B1

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 0.01 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Background

Ephrin-B1 also known as EFNB1, is a member of the ephrin family. The transmembrane-associated ephrin ligands and their Eph family of receptor tyrosine kinases are expressed by cells of the SVZ. Eph/ephrin interactions are implicated in axon guidance, neural crest cell migration, establishment of segmental boundaries, and formation of angiogenic capillary plexi. Eph receptors and ephrins are divided into two subclasses, A and B, based on binding specificities. Ephrin subclasses are further distinguished by their mode of attachment to the plasma membrane: ephrin-A ligands bind EphA receptors and are anchored to the plasma membrane via a glycosylphosphatidylinositol (GPI) linkage, whereas ephrin-B ligands bind EphB receptors and are anchored via a transmembrane domain. An exception is the EphA4 receptor, which binds both subclasses of ephrins.

## Basic Information

### Description

Recombinant Human Ephrin-B1/EFNB1 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Leu 28 - Gly 232 ) of human Ephrin-B1 (Accession #NP\_004420.1) fused with an Fc, 6×His tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Human Ephrin-B1/EFNB1  
Protein was determined by SDS-PAGE with  
Coomassie Blue, showing a band at 65-70  
kDa.