

RP00259

Leader in Biomolecular Solutions for Life Science



# Recombinant Human B7-H5/Gi24/VISTA Protein

Catalog No.: RP00259

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	64115	AAH20568.1

### Tags

C-His

### Synonyms

VSIR;B7-H5;B7H5;C10orf54;DD1alpha;Gi24;PD-1H;PP2135;SISP1;VISTA

## Product Information

Source	Purification
	> 95% by SDS-PAGE.

### Endotoxin

< 0.1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

V-type immunoglobulin domain-containing suppressor of T-cell activation(VISTA) also known as platelet receptor Gi24 , stress-induced secreted protein-1 (Sisp-1) and B7

## Basic Information

### Description

Recombinant Human B7-H5/Gi24/VISTA Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Phe 33 - Ala 194 ) of human VISTA/B7-H5/PD-1H (Accession #AAH20568.1) fused with a 6×His tag at the C-terminus.


### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

## Validation Data

---



Recombinant Human B7-H5/Gi24/VISTA  
Protein was determined by SDS-PAGE with  
Coomassie Blue, showing a band at 30-45  
kDa.