

RP00289

Leader in Biomolecular Solutions for Life Science



Recombinant Human ICAM-2/CD102 Protein

Catalog No.: RP00289

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	3384	P13598

Tags

C-hFc&His

Synonyms

ICAM2;CD102

Product Information

Source	Purification
HEK293 cells	> 97% by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Contact



www.abclonal.com

Background

Basic Information

Description

Recombinant Human ICAM-2/CD102 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Lys25-Gln216) of human ICAM-2 (Accession #NP_000864.2) fused with an Fc, 6xHis tag at the C-terminus.


Bio-Activity

Measured by the ability of the immobilized protein to support the adhesion of THP-1 cells. When cells are added to ICAM-2 coated plates (12.5 µg/mL, 100 µL/well), approximately >40% cells will adhere specifically after 60 minutes at 37 °C.

Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human ICAM-2/CD102 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60-70 kDa.