

RP00410

Leader in Biomolecular Solutions for Life Science



# Recombinant Human IL-5 Protein

Catalog No.: RP00410

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	3567	P05113

### Tags

C-6×His

### Synonyms

IL5;EDF;IL-5;TRF

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

This protein belongs a cytokine that acts as a growth and differentiation factor for both B cells and eosinophils. The encoded cytokine plays a major role in the regulation of eosinophil formation, maturation, recruitment and survival. The increased production of this cytokine may be related to pathogenesis of eosinophil-dependent inflammatory diseases. This cytokine functions by binding to its receptor, which is a heterodimer, whose beta subunit is shared with the receptors for interleukine 3 (IL3) and colony stimulating factor 2 (CSF2/GM-CSF). This gene is located on chromosome 5 within a cytokine gene cluster which includes interleukin 4 (IL4), interleukin 13 (IL13), and CSF2 . This gene, IL4, and IL13 may be regulated coordinately by long-range regulatory elements spread over 120 kilobases on chromosome 5q31.

## Basic Information

### Description

Recombinant Human IL-5 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Ile20-Ser134) of human IL-5 (Accession #P05113) fused with a 6×His tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant protein Human IL-5 was determined by SDS-PAGE under reducing conditions with Coomassie Blue, showing a band at 20 kDa.