

# Recombinant Human TNFRSF11A/CD265 Protein

Catalog No.: RP00458 **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Human	8792	Q9Y6Q6

### Tags

C-His

### Synonyms

TNFRSF11A;CD265;FEO;LOH18CR1;ODF R;OFE;OPTB7;OSTS;PDB2;RANK;TRANC ER

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Background

This protein is a member of the TNF-receptor superfamily. This receptors can interact with various TRAF family proteins, through which this receptor induces the activation of NF-kappa B and MAPK8/JNK. This receptor and its ligand are important regulators of the interaction between T cells and dendritic cells. This receptor is also an essential mediator for osteoclast and lymph node development. Mutations at this locus have been associated with familial expansile osteolysis, autosomal recessive osteopetrosis, and Paget disease of bone. Alternatively spliced transcript variants have been described for this locus.

## Basic Information

### Description

Recombinant Human TNFRSF11A/CD265 Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Ile30-Pro212) of human RANK/TNFRSF11A/CD265 (Accession #Q9Y6Q6) fused with a 6xHis tag at the C-terminus.

### Bio-Activity


### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

## Contact

## Validation Data

---



Recombinant protein Human RANK/TNFRSF11A/CD265 was determined by SDS-PAGE under reducing conditions with Coomassie Blue, showing a band at 30 kDa.