

RP00717

Leader in Biomolecular Solutions for Life Science



Recombinant Mouse B7-H3/CD276 Protein

Catalog No.: RP00717

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Mouse	102657	Q8VE98

Tags

C-His

Synonyms

CD276 antigen; CD276; B7 homolog 3; B7-H3;CD276

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

Background

CD276, also known as B7-H3, is a member of the B7 superfamily with signature IgV and IgG regions in extracellular domains. It is a type I transmembrane protein and shares 20-27% amino acid identity with other B7 family members. B7-H3 is involved in the activation of T lymphocytes, and regulates murine bone formation. It is also reported that B7-H3 may play an important role in muscle-immune interactions, providing further evidence of the active role of muscle cells in local immunoregulatory processes. B7-H3 is expressed on T-cells, natural killer cells, and antigen presenting cells, as well as some non-immune cells, such as osteoblasts, fibroblasts, fibroblast-like synoviocytes and epithelial cells. High expression of B7-H3 in tumor vasculature also correlates with poor survival in patients, suggesting that it may play a role in tumor cell migration.

Basic Information

Description

Recombinant Mouse B7-H3/CD276 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Val29-Phe244) of mouse B7-H3/CD276 (Accession #Q8VE98) fused with a 6×His tag at the C-terminus.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Mouse B7-H3/CD276 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.