

RP00927

Leader in Biomolecular Solutions for Life Science



# Recombinant Human TNFRSF4/OX40/CD134 Protein

Catalog No.: RP00927

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	7293	P43489

### Tags

N-Fc

### Synonyms

TNFRSF4;ACT35;CD134;IMD16;OX40;TXG P1L

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4. Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

## Basic Information

### Description

Recombinant Human TNFRSF4/OX40/CD134 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Leu29-Ala216) of human TNFRSF4/OX40/CD134 (Accession #P43489) fused with an Fc tag at the N-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Human TNFRSF4/OX40/CD134  
Protein was determined by SDS-PAGE under  
reducing conditions with Coomassie Blue.