

RP00946

Leader in Biomolecular Solutions for Life Science



Recombinant Human CTLA-4/CD152 Protein

Catalog No.: RP00946

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	1493	P16410

Tags
C-GST

Synonyms

CTLA4;ALPS5;CD;CD152;CELIAC3;CTLA-4
;GRD4;GSE;IDDM12

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4. Contact us for customized product form or formulation.

Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

Background

Basic Information

Description

Recombinant Human CTLA-4/CD152 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Lys36-Asp161) of human CTLA-4/CD152 (Accession #P16410) fused with a GST tag at the C-terminus.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human CTLA-4/CD152 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.