

RP00964

Leader in Biomolecular Solutions for Life Science



# Recombinant Human IL-1RAcP/IL-1R3 Protein

Catalog No.: RP00964

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	3556	Q9NPH3-2

### Tags

C-Fc & 6×His

### Synonyms

IL1RAP;C3orf13;IL-1RAcP;IL1R3

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.2 μm filtered solution of PBS,pH7.4.Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

## Basic Information

### Description

Recombinant Human IL-1RAcP/IL-1R3 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ser21-Gln356) of human IL-1RAcP/IL-1R3 (Accession #Q9NPH3-2) fused with an Fc, 6×His tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term.<br>After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)